

The Hely & Weber Global Knee™

Addresses Pain, Improves Lives

Surgery is often delayed for Bi & Tricompartmental OA patients:
*Conservative care guidelines, patient logistical needs,
necessary medical and behavioral improvements
and the emotional readiness of the patient often cause delays.*

The Global OA addresses patient pain, provides stability and
assists in mobility, helping patients through these times.



Hybrid Wrap Around Format
Closed calf section with a wrap
around thigh section allows
patients to adjust the size of
the thigh to accommodate
changes in girth.



**Soft, Comfortable
Condyle Pads**
Soft, oversized Comfort
Foam pads accommodate
the inflamed bulbous knee.



**Thin Removable Trakaderm™
"U" Buttress**
This is a "positional buttress"
which can be adjusted for each
individual patient's need. In
addition it can be removed if not
required as a therapeutic modality.

Patella Sling
The "Sling" allows for
directional control of the patella
and offers a slight tilt adjustment.

View from inside of brace reveals extraordinary
craftsmanship and attention to detail.



Our "Patella Sling" controls the
patella with an inferior to superior
directional pull.

**Introducing
our Dynamic
Patella Sling**

*By controlling the
patella in a more
natural position in
the superior femoral
groove, patients can
now manage end range
gait and are able to
propel up and down
stairs with less pain
and more stability.*



**Soft Durable Lycra
Pre-Formed
Popliteal Panel**

*A thin, soft fabric in the
popliteal area provides
more comfort. This results
in less skin irritation, an
overall better fit and
higher patient compliance.*

Functional knee stability
patients haven't seen in years!



Reinforced hinge pockets provide durability and double as a built in pull up assist.

Constructed of breathable lycra urethane laminate and spacer fabric.

The ROM hinges are slightly offset to accommodate inflamed bulbous knees.

Bending and adjusting the struts is made easy with bending irons. This allows hinges to better fit misshapen knees and apply slight varus or valgus force to the knee.



* Bending irons sold separately

Customizable Range of Motion



SIZING CHART

Global Knee™ Part # 7640 PDAC Approved L1832

SIZE	CALF	THIGH
XS	11.0" - 13.0"	15.0" - 17.0"
S	12.5" - 13.5"	16.0" - 18.0"
M	13.5" - 15.5"	17.5" - 19.5"
L	14.5" - 16.5"	18.5" - 20.5"
XL	16.5" - 19.0"	19.5" - 22.0"
XXL	18.0" - 20.5"	20.5" - 23.0"
XXXL	19.0" - 22.0"	24.0" - 27.0"

Larger and custom sizes for 7640 models available upon request. Calf measured at widest part of calf. Thigh measured 6 inches above patella.



See our application video

orders: (800) 654-3241
order@hely-weber.com
www.hely-weber.com



All content © 2015 Hely & Weber

The Hely & Weber
GLOBAL KNEE™

Patent Pending



Designed specifically to address "Global" Knee Osteoarthritis

