

Empire TF Össur Leg Order Package

Send all casts, sockets, files to: 9545 S John Young Parkway, Ste #100, Orlando, FL 32819
tel: (888)839-6213, fax: (800)788-9878 | cs@ossur.com | www.ossur.com

Billing & Shipping Information

Empire's Össur Account # 470109
Contact for ordering or delivery questions:
Contact: Empire Orders
Email: po@empire-medical.com
Phone: 541-245-6657

Ship To:
Clinic:
Address:
City/State/Zip:
CP:
CP Mobile Phone: Fax:
CP email:



Empire PO Number: Shipping Check Priority: Next Day Air **2nd Day**
If no PO number, send quote to Empire at po@empire-medical.com Include return shipping label (\$50 charge inside 48 states)

Completion of this order form with the most accurate and up-to-date information, including all patient information, measurements, and construction information enables us to provide the highest quality product for you and your patient.

Required Patient Information

Last Name: Height: ft. in. Impact Level: High Med Low Amputation Level: TF KD
First Name: Weight: lbs. K-Level: K1 K2 K3 K4 Left Right Bilat Order

Measurements Call to discuss before design

MEASUREMENTS**	SHAPE	DISTAL END SHAPES
 <p>Length: Level Tight* Loose* 5cm 10cm 15cm 20cm 25cm 30cm</p>	<p>Shape by Measurements Shape by cast Shape from Socket Merge scan to Brim Fully Modified CAD file CAD file name:</p> <hr/> <p style="text-align: center;">BRIM</p> <p>Soft IC IC Quad</p> <hr/> <p style="text-align: center;">LINER USED</p>	<p style="text-align: center;">Cylindrical Squared Bulbous Conical</p> <p style="text-align: center;">  </p> <hr/> <p style="text-align: center;">*ALIGNMENT</p> <p style="text-align: center;">Flexion: Adduction:</p> <p>*If a specific alignment isn't provided, ÖCS' standard alignment is 3° flexion & 3° adduction when using expulsion/elevated vacuum suspension or 0° flexion & 0° adduction when using locking suspension.</p> <p>NOTES:</p>
<p>**Length and circumferences are required if Össur Custom Solutions is designing the socket. If not provided and socket fit requires remake, customer is responsible for cost of remake.</p> <p>Cast sent to Össur Custom Solutions (optional) Check here to receive foam carving for this order.</p>	<p>Total Reduction % Ply</p> <p>*The reduction is average of the tight and loose circumference measurements, if both sets are provided.</p>	

Empire TF Össur Leg Order Package

Special Fabrications

Flexible Inner Socket Plastics <i>Suggested L-code: 5651</i>	Fab Custom Inner From Sheet Material <i>Suggested L-codes: 5658 or 5656</i>	Other <i>Suggested L-codes: 5658 or 5656</i>
Flexible EVA with Silicone (Proflex with Silicone)	Lightweight Polyethylene – Firm 25 Shore A (PeLite)	Stovepipe
Flexible EVA without Silicone (Proflex without Silicone)	12mm (1/2") beige white	Slit in lateral side
Flexible EVA – Black (Optek Flex Comfort - Black)	5mm (3/16") beige	Slit in medial side
Flexible EVA – Optek Flex Comfort White + Black (BiLam)	3mm (1/8") beige	
Modified Polyethylene (MPE)		Lightweight Ethylene Vinyl Acetate (Keasy Cone) 6.5mm profile, 10mm distal
EVA (ÖssurFlex)		

Fab Custom Pad

RevoFit

<p>Fab Custom Pad on Check Socket</p> <p>Fab Custom Pad on Laminated Socket</p> <p>Lightweight Polyethylene – Medium Shore A (Plastizote)</p> <p>12mm (1/2") pink white</p> <p>6mm (1/4") pink</p> <p>5mm (3/16") pink</p> <p>Lightweight Polyethylene – Firm 25 Shore A (PeLite)</p> <p>12mm (1/2") beige white</p> <p>5mm (3/16") beige</p> <p>3mm (1/8") beige</p> <p>Location</p> <p>Distal end (1/2") Other (see notes below)</p>	<p>RevoFit on Check Socket <i>Suggested L-code: 5783</i></p> <p>RevoFit on Laminated Socket <i>Suggested L-code: 5783</i></p> <p>Dial Location</p> <p>Distal anterior medial Proximal lateral</p> <p>Other (see notes below)</p> <p>Panel Locations</p> <p>M/L panels (two panels) with Pads</p> <p>Anterior & Posterior panels (two panels) with Pads</p> <p>Anterior only (one panel) with Pad</p> <p>Posterior only (one panel) with Pad</p> <p>Pad on Panel</p> <p>Firm Ethylene Vinyl Acetate 55 Shore A (Puff)</p> <p>12mm (1/2") black</p> <p>10mm (3/8") black</p> <p>6mm (1/4") black</p> <p>Posterior "U" channel with proximal adjustability</p> <p>See markings on cast/socket</p> <p>Other (see notes below)</p>
--	--

Windows and Doors

<p>Windows and Doors on Check Socket</p> <p>Windows and Doors on Laminated Socket</p> <p>Fab Window(s) or cut out(s) Number:</p> <p>Fab panel</p> <p>Fab overlap door (draw window or door on cast/check socket, or use notes field below)</p> <p>Add strap</p> <p>Dacron Velcro</p>	<p>Cosmetics</p> <p>Decorative Lamination on Definitive Using Fred's Legs Sleeve</p> <p>Customer supplied</p> <p>Össur to purchase Fred's item:</p> <p>Customer shaped foam cover (measurements in TF measurement chart) <i>Suggested L-code: 5705 or 5706</i></p> <p>Fab cosmetic distal end</p>
--	--

Notes

Empire TF Össur Leg Order Package

Products

1. Assembly (check one)

Partial Assembly: Customer cuts pylon to final height (default option)



Partial Assembly: Össur cuts pylon to height using landmark to floor measurement



Full Assembly by Customer



Include Össur Leg Bag

2: Liners - Two are ALWAYS Included with Each Leg (check one then fill in desired sizes and options)

Send liners first for measurements/scanning/casting

Suggested L-code: 5673

TF Locking



LINER SIZE

Standard
Conical

Suggested L-code: 5679

Seal-In X5 TF



4Seal Classic



Seal-In TF



LINER SIZE

Standard
Conical

*Suggested L-code: 7700
(only billed with L-5679)*

Seal-In X TF
with Ring






Classic Seal
Grip Seal
Volume Seal

Ring Size:	20	22	24	26	28
	30	32	35	38	41
	47	51	55	60	65

Empire TF Össur Leg Order Package

Products

3: Sockets (check options) No Socket		
Check socket ONLY package	Check socket + Definitive package	Definitive ONLY package
 <p>Polyethylene terephthalate glycol (PETG) High Impact Polystyrene (Orfitrans Stiff, Thermolyn)</p> <p>Flexible Inner Socket Flexible Inner Socket - Plastic Flexible Inner Socket - Sheet Material/Other (PeLite, Keasy Cone) Include second check socket</p>	 <p>Polyethylene terephthalate glycol (PETG) High Impact Polystyrene (Orfitrans Stiff, Thermolyn)</p> <p>Flexible Inner Socket Flexible Inner Socket - Plastic Flexible Inner Socket - Sheet Material/Other (PeLite, Keasy Cone) Include second check socket</p>	 <p>Laminated, with flexible inner Flexible Inner Socket - Plastic Flexible Inner Socket - Sheet Material/Other (PeLite, Keasy Cone) Laminated no flexible inner</p>
<p>NOTE: Exact Material selection is on definitive order form and not needed until the definitive form is ready for submission</p>		

Direct Socket TF

Brim Size:	36	38	40	42	44	46	48.5	51	53.5	56	58.5	61	64
	67	70											
							Fiberglass black resin		Fiberglass neutral resin		Basalt		Carbon
5" Material Kit (standard distal connector, 200 ml resin)													
7" Material Kit (400 ml resin) Medium Distal Large Distal (Most Common)													
9" Material Kit - Large Distal (400 ml resin)													

4. Optional: Connect TF Package Kit

Connect TF Kit

Select your Size:

- Medium Short
- Medium Standard
- Large Short
- Large Standard

Left Right




Note: Connect TF is only available as an option when purchased with Icross TF Locking liner(s), AKOS, and Relax. One Connect TF socket kit per order. Connect TF includes integrated 200 series clutch lock.

Empire TF Össur Leg Order Package

Products


5: Suspension & Socket Adapter - For Conventional Socket Fabrication

Locking Options

Icelock 621 Ratchet & Extra Pin 



Suggested L-code: 5671

Icelock 651 Smooth & Extra Pin 



Suggested L-code: 5671

Icelock 600XM



Suggested L-code: 5671

Icelock 211 Clutch & Extra Pin
Icelock 214 4-Hole Clutch & Extra Pin



Suggested L-code: 5671

Icelock 234 Lanyard



Suggested L-code: 5671

Works With

Icelock 672 Male Pyramid Ti 
Icelock 673 Male Pyramid Al



Icelock 674 4-Hole Adapter Al
Icelock 675 4-Hole Adapter Ti



Icelock Clutch Options for TF

214 AK for left leg with release button on medial side

Suggested L-code: 5671



214 AK for left leg with release button on lateral side

Suggested L-code: 5671



214 AK for right leg with release button on lateral side

Suggested L-code: 5671



214 AK for right leg with release button on medial side

Suggested L-code: 5671



Expulsion Options

Icelock 544 Expulsion Kit for Seal-In



Suggested L-code: 5652

Icelock 552 Expulsion Valve 



Suggested L-code: 5652

Unity Elevated Vacuum Options

Unity Valve TF 



Icelock 544 Unity Plate



Works With

Icelock 544 Socket Adapter

3-Prong Socket Adapter
4-Prong Socket Adapter

Unity TF Tube Kit 



Bulldog Locks

3GEN-A – standard push lock
3GEN-AGG – push lock with grooves
BD5G-63 – ball bearing safety lock
HDT-AP-TL Heavy Duty Top Loading Attachment Plate only (Does not include top load lock, HD pyramid, HD pylon, or HD tube clamp. Those items can be ordered through Bulldog Tools directly or a distributor)

KISS Suspension

KS1-A (Select either a 544 socket adapter or a prong adapter above)
KS2-A
KS4-A

ESP

Lyn Valve - AK RV Auto Expulsion PA0002/B

Empire TF Össur Leg Order Package

Products

6. Direct Socket Suspension Options (Choose when fabricating your socket using a Direct Socket Kit)

Locking

621 Ratchet Pin *Suggested L-code: 5671*



Requires: Modular Spacer Plate



651 Smooth Pin *Suggested L-code: 5671*



Requires: Modular Spacer Plate



214 Clutch Lock *Suggested L-code: 5671*



Add Extra Pin: Long Standard Short

Expulsion

Expulsion Valve 551
Suggested L-code: 5671
Expulsion Valve 552
Suggested L-code: 5671
Lyn Valve - AK RV Auto
Expulsion PA0002/B



Requires: Bottom Seal



Expulsion Plate
Suggested L-code: 5671



Requires: DS Adapter Plate



Unity

Unity TT Valve
Unity TF Valve



Requires: Bottom Seal



544 Unity Plate



Unity TF Tube Kit



Requires: DS Adapter Plate



Waterproof rated Össur products provide permanent protection against harmful ingress of fresh, salt or chlorinated water, even after submersion.

Empire TF Össur Leg Order Package

7. Endoskeletal Components (check all that apply for a single prosthesis)					
4-Hole Adapters		Double Adapters		Single Adapters	
	4-Hole Male		Female Double Adapter 32mm		Male Single Adapter Short
	4-Hole Female		Female Double Adapter 45mm		Male Single Adapter Long
	4-Hole Male w/ Rotation		Female Double Adapter 60mm		Female Single Adapter Short
	4-Hole Male w/ 8mm Ax Offset		Female Double Adapter 75mm		Female Single Adapter Long
	4-Hole Female w/ Rotation	Pylon Kits		Other:	
	4-Hole Female Pyramid 20/20 Offset 10° Right		Female Kit Short	Pylons	
	4-Hole Female Pyramid 20/20 Offset 10° Left		Female Kit Long		Female Short
Prong Adapters		Height Adjustable Pylons		Tube Clamps	
	4-Prong Socket Adapter		Standard		Male Tube Clamp
	3-Prong Socket Adapter		Long		Female Tube Clamp
	Male Pyramid Insert for Prong	Torsion Shock Adapter (Extra Cost)			4-Hole Tube Clamp
	Female Pyramid Insert for Prong		Torsion Shock Adapter <i>Suggested L-codes: 5984 and 5988</i>	<p>* Patients under 220lb will be supplied with Aluminium adaptors unless specified.</p> <p>* Patients between 220-365lb will be provided with Titanium adaptors.</p> <p>Please check this box if you would prefer titanium adaptors.</p>	
	Female Pyramid Insert for Prong LP	Rotation Adapter (Extra Cost)			
			Knee Rotation Adapter (not included in package price) <i>Suggested L-code: 5926</i>		





























Össur may select different components to fit your dimensions if necessary. If you do not want Össur to make any changes, check this box:

8: Knee (check one)		*Knee locking/release cable guides will be laminated into the socket if provided to ÖCS unless requested differently on an order form.			
	<i>Suggested L-codes: 5811</i> Locking Knee		<i>Suggested L-codes: 5611, 5845 and 5925</i> Balance Knee OFM1 w/ IKF Adapter w/ Male Pyramid w/ Loop Adapter		<i>Suggested L-codes: 5812, 5850 and 5925</i> Balance Knee OFM2
	<i>Suggested L-codes: 5814, 5845 and 5850</i> Total Knee 2000 A835300 Pyramid Top Extension Assist option		<i>Suggested L-codes: 5840, 5845 and 5850</i> OHP3 Knee w/ IKF Adapter w/ Male Pyramid w/ Loop Adapter		<i>Suggested L-codes: 5840, 5845 and 5850</i> Paso Knee
	<i>Suggested L-codes: 5828, 5845 and 5925</i> Mauch Knee Mauch Knee Plus†		<i>Suggested L-codes: 5828, 5845, 5848, 5850, 5856 and 5925</i> Rheo Knee		<i>Suggested L-codes: 5828, 5845, 5848, 5850, 5856 and 5925</i> Rheo Knee XC
			<i>Suggested L-codes: 5611, 5845 and 5850</i> Total Knee 1900 A835300 Pyramid Top Extension Assist option		<i>Suggested L-codes: 5814, 5845 and 5850</i> Total Knee 2100 A835300 Pyramid Top Extension Assist option
			<i>Suggested L-codes: 5828, 5845, 5848, 5856 and 5859</i> Power Knee		

Empire TF Össur Leg Order Package

Products

9: Foot: With or Without Unity (check one and fill out the options below)

No Foot		Foot Options:		Foot Size:	cm.	Category (1-9)	Left	Right	Foot Shell:	Beige	Brown		
						Suggested L-Codes†							Suggested L-Codes†
	Balance Foot S  with Unity	5972, 5986	5781*		Balance Foot S Torsion  with Unity	5972, 5984, 5988	5781*						
	Balance J with Unity	5976	5781*		K2 Sensation with DP Flexion with Unity	5972, 5986	5972, 5968					5781*	
	Pro-Flex ST  with Unity 	5981, 5986	5781*		Assure with Unity	5981					5781*		
	Pro-Flex LP Torsion with Unity	5981, 5984, 5988	5781*		Vari-Flex with Unity	5981‡					5781*		
	Pro-Flex Pivot with Unity	5999, See Össur reimbursement guide			Pro-Flex LP  with Unity 	5981, 5986					5781*		
	Pro-Flex LP Align	5981, 5986, 5990			Proprio Foot  with Unity 	5973‡					5781*		
	ReFlex Shock with Unity	5987‡	5781*		ReFlex Rotate with Unity	5987‡, 5984‡					5781*		
	Pro-Flex XC  with Unity 	5981‡	5781*		Pro-Flex XC Torsion with Unity	5987‡, 5984‡					5781*		
	Pro-Flex Modular 	5980‡, 5986‡		Desired build height (mm): With foot cover on. No return or exchange after foot is cut.									
							Foot Size	Build Height					
							22-24cm	Min	Max				
							25-27cm	157mm	313mm				
							28-30cm	168mm	327mm				
								176mm	340mm				

 Waterproof rated Össur products provide permanent protection against harmful ingress of fresh, salt or chlorinated water, even after submersion.

*Unity code added to other codes, e.g. Proprio Foot with Unity: 5973, 5781.

†The responsibility for accurate coding lies solely with the healthcare provider, Össur assumes no responsibility or liability for the provider's coding decisions.

‡PDAC approved.

CLEAR FORM

8 of 10

FOLLOW ÖSSUR ON



ÖSSUR CUSTOM SOLUTIONS
9545 S JOHN YOUNG PARKWAY, STE #100 | ORLANDO, FL 32819
TEL (888) 839-6213, FAX (407) 367-3695
CS@OSSUR.COM | OSSUR.COM



Send all casts, sockets, files to: 9545 S John Young Parkway, Ste #100, Orlando, FL 32819 | tel: (888)839-6213, fax: (800)788-9878 | cs@ossur.com | www.ossur.com

Billing & Shipping Information

Empire's Össur Account # 470109
 Contact for ordering or delivery questions:
 Contact: Empire Orders
 Email: po@empire-medical.com
 Phone: 541-245-6657

Ship To:
 Clinic:
 Address:
 City/State/Zip:
 CP:
 CP Mobile Phone: Fax:
 CP email:


Empire PO Number: If no PO number, send quote to Empire at po@empire-medical.com	Shipping Check Priority: Next Day Air 2nd Day Include return shipping label (\$50 charge inside 48 states)
---	--

Completion of this order form with the most accurate and up-to-date information, including all patient information, measurements, and construction information enables us to provide the highest quality product for you and your patient.

Patient Information

Last Name:	First Name:	Weight:	lbs.	Impact Level: High Medium Low
Amputation Level: AK KD	K-Level: K1 K2 K3 K4	Left: Right: Bilat Order:		

Definitive Socket ■ Call to discuss before fabrication



From Socket From Cast From CAD/AOP file dated:
 Fully Modified CAD file

Lamination

Carbon	Basalt	Heavy Duty Lamination
Decorative Lamination from Fred's Legs Sleeve		Color (PRS Pigment 1-16)
Customer supplied		
Össur to purchase Fred's item:		

Alignment

Transfer	Use Angles
Transfer and zero out	Flexion
Use Lines	Adduction
Other	

Modifications

Socket Sizing				Prosthetic Height (optional)		
Increase	Decrease	Ply:	%:	Transfer	Extend:	Shorten:

Special Fabrications

Flexible Inner Socket Plastics <i>Suggested L-code: 5651</i>	Fab Custom Inner From Sheet Material <i>Suggested L-codes: 5658 or 5656</i>	Other <i>Suggested L-codes: 5658 or 5656</i>
Flexible EVA with Silicone (Proflex with Silicone) Flexible EVA without Silicone (Proflex without Silicone) Flexible EVA – Black (Optek Flex Comfort - Black) Flexible EVA – Optek Flex Comfort White + Black (BiLam) Modified Polyethylene (MPE) EVA (ÖssurFlex)	Lightweight Polyethylene – Firm 25 Shore A (PeLite) 12mm (1/2") beige white 5mm (3/16") beige 3mm (1/8") beige	Stovepipe Slit in lateral side Slit in medial side Lightweight Ethylene Vinyl Acetate (Keasy Cone) 6.5mm profile, 10mm distal

Fab Custom Pad	RevoFit <i>Suggested L-code: 5783</i>
Lightweight Polyethylene – Medium Shore A (Plastizote) 12mm (1/2") pink white 6mm (1/4") pink 5mm (3/16") pink Lightweight Polyethylene – Firm 25 Shore A (PeLite) 12mm (1/2") beige white 5mm (3/16") beige 3mm (1/8") beige Location Distal end (1/2") Other (see notes below)	Dial Location Distal anterior medial Proximal lateral Other (see notes below) Panel Locations M/L panels (two panels) with Pads Anterior & Posterior panels (two panels) with Pads Anterior only (one panel) with Pad Posterior only (one panel) with Pad Pad on Panel Firm Ethylene Vinyl Acetate 55 Shore A (Puff) 12mm (1/2") black 10mm (3/8") black 6mm (1/4") black Posterior "U" channel with proximal adjustability See markings on cast/socket Other (see notes below)

Windows and Doors

Fab Window(s) or cut out(s) Number: Fab panel Fab overlap door (draw window or door on cast/check socket, or use notes field below) Add strap Dacron Velcro	
---	--

Cosmetics

Decorative Lamination on Definitive Using Fred's Legs Sleeve Customer supplied Össur to purchase Fred's item: Customer shaped foam cover (measurements in TF measurement chart) <i>Suggested L-codes: 5705 or 5706</i> Fab cosmetic distal end	
--	--

Notes