

Arctic Hip Brace



CORFLEX®

Designed for life.

CORFLEX.COM

The Arctic Hip brace is universal for ease of use and application while reducing inventory requirements.



The Arctic Hip Brace features a universal and easy-to-fit design with seven points of adjustment for superior individualized fit. Low-profile hinge allows for controlled hip flexion and extension in 15° increments. Two abduction/adduction adjustment points combine with rigid plastic plates to deliver counterforce pressures that help maintain the femoral head in the desired position. Overall brace length is easily adjusted at both the waist and thigh plates, using a secure fastening system that holds its settings for the duration of the treatment protocol.

High-quality soft goods help ensure patient comfort and compliance, featuring padded panel covers and a universal waist belt that accommodates circumferences up to 60". Unique belt closure system also incorporates webbing cinch straps that make it effortless to achieve the desired compression. Easy-release hook closures aid in brace adjustment and removal. Center thigh panel also includes anterior and posterior thumbscrews that allow for targeted offloading of the incision site. A separate sleeve is included with every brace that holds either recirculating cold therapy pads or a 9" x 11" gel pack to encourage valuable cold therapy to help reduce pain and inflammation.

- Lightweight and low-profile hinge provides sturdy support while offering ROM and abduction/adduction control; cold therapy sleeve allows for immediate post-operative cold application to help reduce inflammation.
- Indicated pre or post-op for hip arthroplasty or revision, arthroscopic repairs, dislocation, dysplasia, or instability, post-op proximal hamstring repairs, labral repairs or resections.
- Constructed of aluminum hinge with stainless steel and brass components; plastic cuffs; polyester/foam soft goods with hook and loop closures; Allen tools also included for adjustments.

Suggested Billing Code L-1686

Hip Orthosis, Abduction Control of Hip Joint, Postoperative Hip Abduction Type, Prefabricated, Includes Fitting and Adjustment



Product Number

UNIV

Arctic Hip Brace

89-1500

Accessories

Gel Pack, 9" x 11"

72-7006

Gel Pack, cs/10

72-7006-010

Fits patient heights from 5'2" - 6'6" and up to 60" waist circumference



CORFLEX®
Designed for life.

CORFLEX | 669 EAST INDUSTRIAL PARK DRIVE | MANCHESTER, NH 03109
CORFLEX.COM | TOLL FREE: 800 426 7353 | FAX: 603 623 4111