



HYDRO STEP FOOT

Walk Barefoot with Confidence!

The **HYDRO STEP** foot features an integrated pyramid design, a high-tech keel, and improved traction soles. Its split-toe design accommodates sandals comfortably.

The reinforced treaded sole provides stability and smooth roll over while taking a shower, walking at the beach or around the pool. Users should wear sandals or shoes upon arrival at the water activity and remove them before entering the water. It's important to reapply footwear before leaving.

The **HYDRO STEP** is not intended for walking extended distances bare foot, as this will accelerate the wear and tear on the soling and void the warranty.

Warranty: 1 Year

Features and Benefits

- Designed with zero heel height
- High Tech composite for optimized movement
- Textured sole for improved grip
- Integrated pyramid
- Split toe design



P/N 1129

Product # 1129	
Sizes:	23-28
Weight of Foot Size 27	781g
Build Height of Foot Size 27	97mm
Heel Height:	0mm
Maximum Patient Weight	100kg / 220lb
Warranty	1 year