

UFO-51

Universal AFO

The UFO-51's patent-pending features provide clinicians with fast assembly and unmatched precision.

Flexible Cuff
heat-moldable calf cuff allows for
a wider range of anatomical sizes

Patent Pending Track System
allows for 32 fast, precise height adjustments

Lateral Strut
position uniquely takes advantage of the
footplate variability

Trimmable Footplate
footplate adjustments range from
a US Women's 5 to Men's 14



PART NUMBER	SIZE	FOOT PLATE LENGTH	BRACE HEIGHT
UFOLT	Universal Left	8.75" - 12"	11.5" - 16"
UFORT	Universal Right		